

TEST DATA SHEET FOR ELECTROSURGICAL UNIT IDS SERIES

W/ O NUMBER: _____

Catalog #: _____ Unit Serial #: _____ Assy # 18-077- & Rev. #: _____

Test using QA-256

6.1.3.1 Final Assembly Router completed up to Step 3

☐ Completed

6.1.3.2 Verify Serial Number on back of Unit to DHR and Final Assembly Router

☐ Completed

6.1.3.3 Check hardware (8 screws on Main PCB, 4 Nuts on Display PCB) are present,
Harnesses are plugged in correctly and Relay PCB is seated securely.

☐ Completed

6.1.3.4 Verify on Main PCB that C40 is not in contact with R42

..... ☐ Completed

6.1.3.5 Verify on Main PCB that U22 is not in contact with legs of TR3

..... ☐ Completed

6.2 NEM TESTING

| | | | |
|---|---|--|-----------------------|
| Pre Burn-in Testing Performed By/Date: _____ | | | |
| NEM Calibrated | <input type="checkbox"/> Pass <input type="checkbox"/> Fail | Split pad upper trip res(6.2.13) $130 \leq \quad \leq 140 \Omega$ | |
| In Solid Pad @ 0Ω | <input type="checkbox"/> Pass <input type="checkbox"/> Fail | Solid pad upper trip res. $5 \leq \quad \leq 9\Omega$ | |
| Solid pad limit once established | $20 \leq \quad \leq 25\Omega$ | Split pad lower trip res. $8 \leq \quad \leq 12\Omega$ | |
| Split pad low limit | $2 \leq \quad \leq 6\Omega$ | TEST FIXTURES | |
| % increase in split to cause alarm | _____ % | Resistance Substitution | Eq. # _____ |
| Upper hard limit res. | $145 \leq \quad \leq 155\Omega$ | Dial Resistance Test Fixture | # <u>F-273-</u> _____ |

6.3 Service Mode Data / 5 VDC

| MULTIMETER ID #: _____ | | | + 5 VDC Reading taken at Min LED setting | | 4.90 $\leq \quad \leq$ 5.10 VDC |
|--------------------------|---------------|-----------|--|------------------------------|---------------------------------|
| Mode | Power Setting | V Setting | V (No Activation) | V (Open Circuit) | V (Arcing) |
| Cut I (Cut) | MAX POWER | -12 | $10.4 \leq \quad \leq 13.6V$ | $10.4 \leq \quad \leq \quad$ | $10.4 \leq \quad \leq 13.6V$ |
| | | 5.0 | $4.2 \leq \quad \leq 6.0V$ | $4.2 \leq \quad \leq 6.0V$ | $4.2 \leq \quad \leq 6.0V$ |
| | | 2.5 | $2.2 \leq \quad \leq 3.3V$ | $2.2 \leq \quad \leq 3.3V$ | $2.2 \leq \quad \leq 3.3V$ |
| | | 3.3 | $3.0 \leq \quad \leq 4.4V$ | $3.0 \leq \quad \leq 4.4V$ | $3.0 \leq \quad \leq 4.4V$ |
| | | 12 | $10.4 \leq \quad \leq \quad$ | $10.4 \leq \quad \leq \quad$ | $10.4 \leq \quad \leq 13.6V$ |
| Spray | 80W | -12 | $10.4 \leq \quad \leq \quad$ | $10.4 \leq \quad \leq \quad$ | $10.4 \leq \quad \leq 13.6V$ |
| | | 5.0 | $4.2 \leq \quad \leq 6.0V$ | $4.2 \leq \quad \leq 6.0V$ | $4.2 \leq \quad \leq 6.0V$ |
| | | 2.5 | $2.2 \leq \quad \leq 3.3V$ | $2.2 \leq \quad \leq 3.3V$ | $2.2 \leq \quad \leq 3.3V$ |
| | | 3.3 | $3.0 \leq \quad \leq 4.4V$ | $3.0 \leq \quad \leq 4.4V$ | $3.0 \leq \quad \leq 4.4V$ |
| | | 12 | $10.4 \leq \quad \leq \quad$ | $10.4 \leq \quad \leq \quad$ | $10.4 \leq \quad \leq 13.6V$ |
| JP8 Jumper Removed | | | (✓ when done) _____ | | |
| Performed By/Date: _____ | | | | | |

| | | |
|----------------------------|---|--|
| 6.4 Continuity Test | | SAFETY/ECG ANALYZER ID #: _____ |
| AC GND to Chassis Ground | $0.10 \geq \quad \Omega$ | |
| 6.5 Line Current | | |
| Cut Line Current | (400W) $8.00 > \quad A$ (others) $4.50 > \quad A$ | |
| Coag Line Current | $1.85 > \quad A$ | |

TEST DATA SHEET FOR ELECTROSURGICAL UNIT IDS SERIES

Unit Serial #: _____

| 6.6 Low Freq. Risk Current | | Normal | Reverse |
|----------------------------|-----------|---|---|
| Chassis Grounding Lug | 500uA MAX | $500 \geq \underline{\hspace{1cm}} \mu\text{A}$ | $500 \geq \underline{\hspace{1cm}} \mu\text{A}$ |
| Return Electrode (1) | 10uA MAX | $10 \geq \underline{\hspace{1cm}} \mu\text{A}$ | $10 \geq \underline{\hspace{1cm}} \mu\text{A}$ |
| Return Electrode (2) | 10uA MAX | $10 \geq \underline{\hspace{1cm}} \mu\text{A}$ | $10 \geq \underline{\hspace{1cm}} \mu\text{A}$ |
| Monopolar Ft-Control | 10uA MAX | $10 \geq \underline{\hspace{1cm}} \mu\text{A}$ | $10 \geq \underline{\hspace{1cm}} \mu\text{A}$ |
| Monopolar Handpiece | 10uA MAX | $10 \geq \underline{\hspace{1cm}} \mu\text{A}$ | $10 \geq \underline{\hspace{1cm}} \mu\text{A}$ |
| Monopolar Handpiece (2) | 10uA MAX | $10 \geq \underline{\hspace{1cm}} \mu\text{A}$ | $10 \geq \underline{\hspace{1cm}} \mu\text{A}$ |
| Monopolar Handpiece (3) | 10uA MAX | $10 \geq \underline{\hspace{1cm}} \mu\text{A}$ | $10 \geq \underline{\hspace{1cm}} \mu\text{A}$ |
| Bipolar Active | 10uA MAX | $10 \geq \underline{\hspace{1cm}} \mu\text{A}$ | $10 \geq \underline{\hspace{1cm}} \mu\text{A}$ |
| Bipolar Return | 10uA MAX | $10 \geq \underline{\hspace{1cm}} \mu\text{A}$ | $10 \geq \underline{\hspace{1cm}} \mu\text{A}$ |
| Mono Footswitch (Pin 1) | 500uA MAX | $500 \geq \underline{\hspace{1cm}} \mu\text{A}$ | $500 \geq \underline{\hspace{1cm}} \mu\text{A}$ |
| Mono Footswitch (Pin 2) | 500uA MAX | $500 \geq \underline{\hspace{1cm}} \mu\text{A}$ | $500 \geq \underline{\hspace{1cm}} \mu\text{A}$ |
| Mono Footswitch (Pin 3) | 500uA MAX | $500 \geq \underline{\hspace{1cm}} \mu\text{A}$ | $500 \geq \underline{\hspace{1cm}} \mu\text{A}$ |
| Mono Footswitch (Pin 4) | 500uA MAX | $500 \geq \underline{\hspace{1cm}} \mu\text{A}$ | $500 \geq \underline{\hspace{1cm}} \mu\text{A}$ |
| Bipolar Footswitch (Pin 1) | 500uA MAX | $500 \geq \underline{\hspace{1cm}} \mu\text{A}$ | $500 \geq \underline{\hspace{1cm}} \mu\text{A}$ |
| Bipolar Footswitch (Pin 2) | 500uA MAX | $500 \geq \underline{\hspace{1cm}} \mu\text{A}$ | $500 \geq \underline{\hspace{1cm}} \mu\text{A}$ |
| Bipolar Footswitch (Pin 3) | 500uA MAX | $500 \geq \underline{\hspace{1cm}} \mu\text{A}$ | $500 \geq \underline{\hspace{1cm}} \mu\text{A}$ |

Performed By/Date: _____ ☐ Pass ☐ Fail

| 6.7 INITIAL OUTPUT POWER TEST | | ES ANALYZER ID # _____ | Output [Watts] |
|-------------------------------|----------------------------------|---|----------------|
| Cut I (Max Power) | 400W @ 500Ω, 300W or 200W @ 300Ω | (400W) $320 \leq \underline{\hspace{1cm}} \leq 480\text{W}$ or (300W) $240 \leq \underline{\hspace{1cm}} \leq 360\text{W}$ or (200W) $160 \leq \underline{\hspace{1cm}} \leq 240\text{W}$ | |
| Cut II (if applicable) | 300W @ 300Ω | $240 \leq \underline{\hspace{1cm}} \leq 360\text{W}$ | |
| Blend | 200W @ 300Ω | $160 \leq \underline{\hspace{1cm}} \leq 240\text{W}$ | |
| Pinpoint | 120W @ 500Ω | $96 \leq \underline{\hspace{1cm}} \leq 144\text{W}$ | |
| Spray | 80W @ 500Ω | $64 \leq \underline{\hspace{1cm}} \leq 96\text{W}$ | |
| Bipolar (Max Power) | 120W or 80W @ 150Ω | (120W) $96 \leq \underline{\hspace{1cm}} \leq 144\text{W}$ or (80W) $64 \leq \underline{\hspace{1cm}} \leq 96\text{W}$ | |

Performed By/Date: _____ ☐ Pass ☐ Fail

6.8 Burn-In: Start Time: ____:____ Stop Time: ____:____ ☐ Pass ☐ Fail

Jumper on JP6 removed (✓ when done): _____

6.8.6 Post Burn-in drift check: Split pad upper trip res. after burn-in _____Ω

Split Change in Resistance (Drift) difference from 6.2.13 and 6.8.6 $0 \leq \underline{\hspace{1cm}} \leq 10\Omega$ ☐ Pass ☐ Fail
Resistance Substitution Box Eq. # _____

| 6.9 HIPOT TEST | | Test Voltage | HIPOT ID # _____ | |
|---|--|--------------|-------------------------------|--|
| AC Line to: Chassis Ground | | 1.84kVAC | <input type="checkbox"/> Pass | <input type="checkbox"/> Fail |
| AC Line to: Return Electrode (Pole 1) | | 3.0kVAC | <input type="checkbox"/> Pass | <input type="checkbox"/> Fail |
| AC Line to: Return Electrode (Pole 2) | | 3.0kVAC | <input type="checkbox"/> Pass | <input type="checkbox"/> Fail |
| AC Line to: Monopolar Foot-Control Active | | 3.0kVAC | <input type="checkbox"/> Pass | <input type="checkbox"/> Fail <input type="checkbox"/> N/A |
| AC Line to: Monopolar Handpiece Active | | 3.0kVAC | <input type="checkbox"/> Pass | <input type="checkbox"/> Fail |
| AC Line to: Bipolar Active | | 3.0kVAC | <input type="checkbox"/> Pass | <input type="checkbox"/> Fail |
| Performed By / Date: _____ | | | | |

TEST DATA SHEET FOR ELECTROSURGICAL UNIT IDS SERIES

Unit Serial #: _____

6.10 High Frequency Leakage Testing, Long Leads

FM Power Meter ID #: _____

| Test @ 200Ω | Active (Foot) | Active (Hand) | Return (Hand) | Return (Foot) | Pass/ Fail |
|-------------|---------------|---------------|---------------|---------------|--|
| Cut I | 150≥_____mA | 150≥_____mA | 150≥_____mA | 150≥_____mA | <input type="checkbox"/> Pass <input type="checkbox"/> Fail |
| Cut II | 150≥_____mA | 150≥_____mA | 150≥_____mA | 150≥_____mA | <input type="checkbox"/> Pass <input type="checkbox"/> Fail <input type="checkbox"/> N/A |
| Blend (Max) | 150≥_____mA | 150≥_____mA | 150≥_____mA | 150≥_____mA | <input type="checkbox"/> Pass <input type="checkbox"/> Fail |
| Pinpoint | 150≥_____mA | 150≥_____mA | 150≥_____mA | 150≥_____mA | <input type="checkbox"/> Pass <input type="checkbox"/> Fail |
| Spray | 150≥_____mA | 150≥_____mA | 150≥_____mA | 150≥_____mA | <input type="checkbox"/> Pass <input type="checkbox"/> Fail |
| Bipolar | 63≥_____mA | N/A | N/A | 63≥_____mA | <input type="checkbox"/> Pass <input type="checkbox"/> Fail |

Performed By/Date: _____

Display, Audio, Errors, Arcing, Cross couple and Activation Errors

| | | | |
|----------------------------|---|--------------------------------|--|
| 6.11 Displays, LED's | <input type="checkbox"/> Pass <input type="checkbox"/> Fail | 6.14 ARC Test (All Modes) | <input type="checkbox"/> Pass <input type="checkbox"/> Fail |
| 6.12 Audio | <input type="checkbox"/> Pass <input type="checkbox"/> Fail | 6.14.11 Cross Coupling Test | <input type="checkbox"/> Pass <input type="checkbox"/> Fail <input type="checkbox"/> N/A |
| 6.13 Accessory Inlet | <input type="checkbox"/> Pass <input type="checkbox"/> Fail | 6.15 Activation Errors F1 – F6 | <input type="checkbox"/> Pass <input type="checkbox"/> Fail |
| Test Fixture #: _____ | Performed By/Date: _____ | | |

6.16 Load Curve

ANALYZER ID#: _____

| Mode | Dial Setting | Pass/Fail | Mode | Dial Setting | Pass /Fail |
|-------------|-----------------|---|----------|--------------|--|
| Cut I | Max W | <input type="checkbox"/> Pass <input type="checkbox"/> Fail | Cut II | 300W | <input type="checkbox"/> Pass <input type="checkbox"/> Fail <input type="checkbox"/> N/A |
| Cut I | 150W | <input type="checkbox"/> Pass <input type="checkbox"/> Fail | Cut II | 150W | <input type="checkbox"/> Pass <input type="checkbox"/> Fail <input type="checkbox"/> N/A |
| Cut I | 50W | <input type="checkbox"/> Pass <input type="checkbox"/> Fail | Cut II | 50W | <input type="checkbox"/> Pass <input type="checkbox"/> Fail <input type="checkbox"/> N/A |
| Blend (Max) | 200W | <input type="checkbox"/> Pass <input type="checkbox"/> Fail | Pinpoint | 120W | <input type="checkbox"/> Pass <input type="checkbox"/> Fail |
| Blend (Max) | 100W | <input type="checkbox"/> Pass <input type="checkbox"/> Fail | Pinpoint | 60W | <input type="checkbox"/> Pass <input type="checkbox"/> Fail |
| Blend (Max) | 50W | <input type="checkbox"/> Pass <input type="checkbox"/> Fail | | | |
| Bipolar | Max W | <input type="checkbox"/> Pass <input type="checkbox"/> Fail | Spray | 80W | <input type="checkbox"/> Pass <input type="checkbox"/> Fail |
| Bipolar | 60W (400W only) | | Spray | 40W | <input type="checkbox"/> Pass <input type="checkbox"/> Fail |
| Bipolar | 40W (200W/300W) | <input type="checkbox"/> Pass <input type="checkbox"/> Fail | | | |

Performed By / Date: _____

6.17 Waveforms

OSCILLOSCOPE ID # _____

HV PROBE ID # _____

| Mode / Setting | Observed waveform matches: | Mode / Setting | Observed waveform matches: |
|------------------|---|----------------|--|
| Cut I @ Max W | 200V/Div 2uSec <input type="checkbox"/> Pass <input type="checkbox"/> Fail | Pinpoint @ Max | 500V/Div 10uSec <input type="checkbox"/> Pass <input type="checkbox"/> Fail |
| Cut II @ Max | 200V/Div 2uSec <input type="checkbox"/> Pass <input type="checkbox"/> Fail <input type="checkbox"/> N/A | Spray @ Max | 1000V/Div 10uSec <input type="checkbox"/> Pass <input type="checkbox"/> Fail |
| Blend Min@ 200W | 500V/Div 10uSec <input type="checkbox"/> Pass <input type="checkbox"/> Fail | Bipolar @ Max | 200V/Div 5uSec <input type="checkbox"/> Pass <input type="checkbox"/> Fail |
| Blend Max @ 200W | 500V/Div 10uSec <input type="checkbox"/> Pass <input type="checkbox"/> Fail | | |

Performed By/Date: _____

TEST DATA SHEET FOR ELECTROSURGICAL UNIT IDS SERIES

Unit Serial #: _____

6.18 Output Power

ES ANALYZER ID # _____

| Mode / Setting | Output Range | Output Power | Mode / Setting | Output Range | Output Power |
|--|--------------|---|----------------|--------------|---------------------|
| CUT I @ 300 Ω (500 Ω for 400W models) | Watts | Watts | BLEND @ 300 Ω | Watts | Watts |
| 6 | 1 - 11 | 1 ≤ _____ ≤ 11 W | 6 | 1 - 11 | 1 ≤ _____ ≤ 11 W |
| 50 | 40 - 60 | 40 ≤ _____ ≤ 60 W | 50 | 40 - 60 | 40 ≤ _____ ≤ 60 W |
| 100 | 80 - 120 | 80 ≤ _____ ≤ 120 W | 100 | 80 - 120 | 80 ≤ _____ ≤ 120 W |
| 150 | 120 - 180 | 120 ≤ _____ ≤ 180 W | 150 | 120 - 180 | 120 ≤ _____ ≤ 180 W |
| 200 | 160 - 240 | 160 ≤ _____ ≤ 240 W | 200 | 160 - 240 | 160 ≤ _____ ≤ 240 W |
| 250 (if applicable) | 200 - 300 | 200 ≤ _____ ≤ 300 W <input type="checkbox"/> | | | |
| 300 (if applicable) | 240 - 360 | 240 ≤ _____ ≤ 360 W <input type="checkbox"/> | | | |
| 400 (if applicable) | 320 - 480 | 320 ≤ _____ ≤ 480 W <input type="checkbox"/> | | | |
| Mode / Setting | Output Range | Output Power | Mode / Setting | Output Range | Output Power |
| Pinpoint @ 500 Ω | Watts | Watts | Spray @ 500 Ω | Watts | Watts |
| 6 | 1 - 11 | 1 ≤ _____ ≤ 11 W | 6 | 1 - 11 | 1 ≤ _____ ≤ 11 W |
| 60 | 48 - 72 | 48 ≤ _____ ≤ 72 W | 40 | 32 - 48 | 32 ≤ _____ ≤ 48 W |
| 120 | 96 - 144 | 96 ≤ _____ ≤ 144 W | 80 | 64 - 96 | 64 ≤ _____ ≤ 96 W |
| Mode / Setting | Output Range | Output Power | | | |
| Bipolar @ 150 Ω | Watts | Watts | | | |
| 6 | 1 - 11 | 1 ≤ _____ ≤ 11 W | | | |
| 40 | 32 - 48 | 32 ≤ _____ ≤ 48 W | | | |
| 80 | 64 - 96 | 64 ≤ _____ ≤ 96 W | | | |
| 120 (if applicable) | 96 - 144 | 96 ≤ _____ ≤ 144 W <input type="checkbox"/> N/A | | | |

Performed By/Date: _____ ☐ Pass ☐ Fail

6.19 Presets ☐ Pass ☐ Fail ☐ N/A

6.19.6 Factory Presets set ☐ Pass ☐ Fail ☐ N/A

6.20 Pre-Close Verification

| | | | | | | |
|---------------------|-------------------------------|---------------------|---------------------------|-------------------------|-----------------------------|----------------------|
| All Jumpers Removed | All Cables Securely Connected | No Foreign Material | Cover on, Invert, & Shake | No Automatic Activation | Set Power Levels to Minimum | Check DHR and Router |
| _____ | _____ | _____ | _____ | _____ | _____ | _____ |

6.21 Final Verification/ Certification and

InfinityQS Data Verification Performed By/ Date: _____

Accessory Kit Part# _____ Accessory Kit Lot# _____

QA Inspection per QA-196 by/ Date: _____

Packaging Completed per WI-0387 By/ Date: _____